

Patient gown - dimensions and tolerances given in [cm]

Parameter		Size					Tolerance
		Universal Plus	Universal	132/146	104/128	54/98	
A	Length of the gown (from the highest shoulder point)	120	110	90	70	50	± 2
B	Gown width (at the bottom)	170	140	120	112	92	± 2
C	Arm length	29	27	18	16	14	± 0,5
D	Length of the ties to tie at the neck	2 x 20	2 x 20	2 x 20	2 x 20	2 x 20	± 1
E	Length of the ties to tie at the waist	210	180	160	140	120	± 5

